

# Luminex Xponent: Plate-level Quality Control

Sahal Thahir & Jeff Bailey

2024-06-04

## Contents

Aim . . . . .	1
---------------	---

## Aim

This code aims to perform plate-level quality control analysis for Luminex studies through bead counts, Background MFI, and Standard curve analysis.

## Instructions

### 1. Modify YAML Parameters:

- Open this document in RStudio.
- Edit the YAML parameters at the top of the document:
  - `input_file`: Path to your input CSV file.
  - `min_beadcount`: Set the minimum bead count threshold.
  - `min_Rsquared`: Set the minimum R-squared value for standard curve QC.
  - `file_path`: Specify the output directory for plots and PDFs.
  - ‘

### 2. Run the Analysis:

- After modifying the YAML parameters, click the *knit* button at the top of the page to generate the report.

## Code Explanation

### Dependencies

- R version 4.3.1 (2023-06-16)
  - `tidyverse` (v. 2.0.0)
  - `here` (v. 1.0.1)
  - `gridExtra` (v. 2.3)

This code chunk will ensure these packages are installed and loaded.

### Import and Data Extraction

#### 1. Raw data import: Import raw data CSV file into R.

```
# No input required
#YAML extraction
# Read parameters from YAML header
input_file <- params$input_file
min_beadcount <- params$min_beadcount
min_Rsquared <- params$min_Rsquared
file_path <- params$file_path
```

```
background_col <- params$background

# Raw data import
lines <- readLines(input_file)
```


```
## Warning in readLines(input_file): incomplete final line found on
## 'data/raw/luminex/Control Run 13 June 2024_20240613_161509.csv'
```

- 2. **Extract dataframes:** find median MFI and bead count data. This includes a check to ensure the median MFI and beadcount dataframes align.
- 3. **Identify low bead count wells:** this searches the beadcount data frame for well-antigen data with low bead count, as set in the YAML.This will export a csv file of the samples which will need to be rerun.

```
## The bead QC list has been exported to data/qc/Control Run 13 June 2024_20240613_161509_beadqc_low_df
```

Visualization of bead count data


This code generates a well schematic identifying wells identified as having low bead counts. as well as a heatmap of antigen/well locations and the quantity of beads, to assist in determining the cause of low bead counts (low antigen concentration in mastermix, vortexing between wells, etc. )



```
## # A tibble: 772 x 4
##   Location Sample      Antigen      Plate
##   <chr>      <chr>      <chr>      <chr>
## 1 1(1,A1) Background0 BSA      Control Run 13 June 2024_20240613_161509
## 2 1(1,A1) Background0 Rickettsia Control Run 13 June 2024_20240613_161509
## 3 1(1,A1) Background0 Pf.MSP1    Control Run 13 June 2024_20240613_161509
```

```
## 4 1(1,A1) Background0 Dengue.3 Control Run 13 June 2024_20240613_161509
## 5 1(1,A1) Background0 Dengue.4 Control Run 13 June 2024_20240613_161509
## 6 1(1,A1) Background0 CCHFV Control Run 13 June 2024_20240613_161509
## 7 1(1,A1) Background0 Dengue.2 Control Run 13 June 2024_20240613_161509
## 8 1(1,A1) Background0 Measles Control Run 13 June 2024_20240613_161509
## 9 1(1,A1) Background0 Chikungunya Control Run 13 June 2024_20240613_161509
## 10 1(1,A1) Background0 Pv.MSP1 Control Run 13 June 2024_20240613_161509
## # i 762 more rows

## The bead QC list has been exported to data/qc/Control Run 13 June 2024_20240613_161509_beadqc_low_df
## The bead QC report has been exported to data/qc/Control Run 13 June 2024_20240613_161509_beadqc.png
```



4. **Modify Median MFI data:** All plate locations with a low beadcount will be switched to a NA value for standard curve evaluation

```
## Median MFI data has been modified for low bead count wells.
```

5. **Evaluate standard curves:** Make sure the background and standard analytes are defined in the YAML params

This first chunk conducts background MFI subtraction, based on the analyte defined in the YAML.

```
## Background MFI subtraction completed.
```

This code chunk

```
## No valid plots to display.
```

```
## [1] Analyte R_squared Equation
```

```
## <0 rows> (or 0-length row.names)
```

### Create quality controlled Median MFI data

```
## [1] Analyte    R_squared Equation  min_Rsq
## <0 rows> (or 0-length row.names)

## [1] "This plate has failed plate-level QC for standard curves."
## [1] "No analytes met the preset minimum R-squared (see min_Rsquared in YAML)."
```

## [1] "All MFI values from this plate will display as 'Failed Plate QC'."

## Plate failed QC. All median MFI data set to 'Failed QC (standards)'.

## The cleaned median MFI data has been exported to data/clean/luminex/Control Run 13 June 2024\_2024061

### Developer notes

->